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Servona Silikonkurzkanüle

Servona Stoma-Button

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Servona Short Silicone Tracheostomy Tube | Servona Stoma-Button

1. General Description

The Servona Short Silicone Tracheostomy Tube and the Servona Stoma-Button are inserted in an surgically created tracheostoma and thus guarantee that the tracheostoma remains stable and open.

The Servona Short Silicone Tracheostomy Tube and the Servona Stoma-Button are high quality products made of medical-grade silicone and are supplied in non-sterile packaging.

The Servona Short Silicone Tracheostomy Tube and the Servona Stoma-Button are not made with natural rubber latex as well as PVC and do not contain any metallic materials, so that there are no restrictions in the use of magnetic resonance imaging (MRI).

The material properties and the good workmanship mean that the Servona Short Silicone Tracheostomy Tube and the Servona Stoma-Button are well tolerated tracheostomy products that can also be used on a long-term basis.

The soft silicone adapts itself exactly to the anatomical conditions in the trachea, thus guaranteeing a high level of wearing comfort.

The Servona Short Silicone Tracheostomy Tube and the Servona Stoma-Button are one-patient products.

2. Indications

The Servona Short Silicone Tracheostomy Tube and the Servona Stoma-Button have been developed for the rehabilitation of patients after a tracheotomy or a laryngectomy. They are inserted into the tracheostoma and guarantee that the tracheostoma remains stable and open.

3. Contraindications

It is advised not to use the Servona Short Silicone Tracheostomy Tube and Stoma-Button in case of known allergies to silicone or an excessively narrow or unsuitably positioned tracheostoma.

4. Side Effects and Warnings

The following side effects may occur when the Servona Short Silicone Tracheostomy Tube or the Servona Stoma-Button is used:

- Irritation of tissue in the region of the tracheostoma through excessive pressure on the surrounding tissue
- Contamination of the surrounding tissue by yeast fungi (candida), bacteria, etc. (see section on Cleaning and Care)
- In the fenestrated version of the Servona Short Silicone Tracheostomy Tube
 - Possible formation of granulation tissue
 - Possible obstruction of the fenestration in the event of intense secretion formation or viscous consistency of the secretion

In case of any questions or complications please consult your physician.

5. Application

Caution: Please do not use oily substances or silicone spray as lubricants, as these may change the material properties of the silicone.

Before inserting the Servona Short Silicone Tracheostomy Tube or the Servona Stoma-Button, we recommend that they are cleaned with a suitable cleaning solution.

Before inserting the device, ensure that you have selected the correct size. When removing the product from the packaging, it is also important to ensure that the packaging is unopened and intact. Before inserting the product, check it for possible damage.

a. Inserting and Removing the Servona Short Silicone Tracheostomy Tube and the Servona Stoma-Button

Servona Short Silicone Tracheostomy Tube:

Before inserting the Servona Short Silicone Tracheostomy Tube in the tracheostoma, ensure that the tracheostomy tube is properly lubricated.

Secure the Servona Short Silicone Tracheostomy Tube with a tracheostomy tube holder. The tracheostomy tube holder must be positioned in such a way that the freedom of movement of the head is not restricted and that a firm hold of the tracheostomy tube in the tracheostoma is ensured.

Please use only recommended tracheostomy tube holders. Tracheostomy tube holders with sharp-edge fastening may damage the retaining eyelets.

Servona Stoma-Button:

To insert the Servona Stoma-Button, press the rear flange of the Servona Stoma-Button together. In this state, the Servona Stoma-Button is inserted in the tracheostoma. The flange develops in the tracheostoma and keeps the Servona Stoma-Button in position.

Removal:

The Servona Short Silicone Tracheostomy Tube and the Servona Stoma-Button are removed by pulling them carefully out of the tracheostoma.

b. Cleaning and Care**Notes:**

The products must not be resterilized and must not be exposed to direct sunlight as this may lead to the destruction of the material.

- To prolong the service life of the Servona Short Silicone Tracheostomy Tube and the Servona Stoma-Button, they should be cleaned at least once a day or more frequently if necessary. Before cleaning, the product must be removed from the tracheostoma.

The Servona Short Silicone Tracheostomy Tube and the Servona Stoma-Button can be heated up to 65° C in a water bath. The

Servona Short Silicone Tracheostomy Tube and the Servona Stoma-Button can be cleaned using mild household cleaning agents. We recommend the use of a cleaning solution for tracheostomy tubes. Rinse the products thoroughly under running water after being cleaned and dry them with a lint-free cloth.

Caution: Under no circumstances may cleaning be carried out with aggressive household detergents or high-percentage alcohol as these may alter the material properties. Furthermore, do not use any disinfectants that release chlorine or agents that contain strong lysers.

Secretions sticking to the tube can be carefully removed using a cleaning brush. Make sure you do not apply excessive pressure. To prevent damaging the silicone, mechanical or sharp-edged aids such as scrapers, knives, fingernails or the like must not be used under any circumstances.

Warning: The cleaning brush must not be inserted into the tracheostoma, as it can be aspirated and lead to suffocation!

After cleaning, the aid must be stored in a dry place.

Even after being cleaned several times, in accordance with fixed standards the Servona Short Silicone Tracheostomy Tube and the Servona Stoma-Button can be used several times for a patient for at least 29 days and a maximum of 6 months, when following all cleaning application instruction.

6. Warnings:

Do not modify your product in any way!

Any adaptations like additional fenestration may be carried out by qualified personnel only!

If the product presents any damage such as sharp edges or cracks, do not use it and proceed to return to your provider.

Non-observance of these instructions can have serious effects on the health of the patient and cause life-threatening complications. The manufacturer does not accept any liability for any damage resulting from non-observance of these instructions or from a manipulation of the product.

7. Safety

Each product is assigned a batch number, which is marked on the outer packaging and on the label, to allow unequivocal traceability of the product. The batch number must be entered, or the enclosed label affixed, in the tracheostomy tube card supplied. It is recommended that the patient should always carry his tracheostomy tube card with him.

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