Tracheostomapolster  
Tracheostomy Stoma Pad  
Coscinetto tracheostomica

DE

EN

IT

# SERVOX<sup>®</sup> Seal Pad

CE 0123 MD

**T&L Co., Ltd.**

70-17, Wonam-ro, Won-gok-myeon  
Anseong-si, Gyeonggi-do  
17554 Korea  
Tel. +82 31 651 5220  
Fax. +82 31 651 6212

**KTR Europe GmbH**

Mergenthalerallee 77,  
65760 Eschborn, Germany

**Orbisana Healthcare GmbH**

Biberweg 2426  
53842 Troisdorf | Germany  
Tel. +49 (0)2241 9322-0  
Fax +49 (0)2241 9322-277  
info@orbisana.de  
www.orbisana.de

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## SERVOX® Seal Pad

### Tracheostomapolster

## 1. Vorwort

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### Tracheostomy Stoma Pad

## 1 Preface

These instructions for use are intended for patients, relatives and medical professionals in the home and clinical setting.

### Note

Read these instructions for use carefully. They provide a lot of important information for safe and proper use. Keep these instructions for future reference.

For traceability, SERVOX® Seal Pads are assigned a batch number. The number is indicated on the packaging following the 'LOT' marking.

### 1.1 Intended use

The SERVOX® Seal Pad absorbs tracheal secretions to reduce irritation caused by heavy secretions at the tracheostoma and cushions the cannula shield to prevent pressure points. It evens out irregularities on the tracheostoma and thus seals leaks.

### 1.2 Indication

Patients using a tracheostomy tube with an outer diameter of approx. 10 to 15 mm.

Patients with uneven tracheostomes tracheostomas.

Patients with a speaking valve or a tracheostoma valve.

### 1.3 Contraindication

Allergies against polyurethane.

### 1.4 Side effects

If the indication, contraindication, safety- and application instructions are observed, no side effects are reported.

## 2 Safety instructions

Not following these safety instructions can result in serious health effects or even life-threatening situations. The manufacturer accepts no liability in the event of damage caused by not following the instructions or by manipulation of the product.

## Risk of infection

- Do not reuse!
- Do not use if sterile packaging is damaged!
- Do not use after expiry date!
- Replace after 24 hours at the latest or more frequently, if necessary!
- In case of heavy secretion, the compress may swell, or the front side may become discoloured. Then if that happens, you must replace the compress!

### Note

We ask you to inform us and, if necessary, the competent authority if you suffer a life-threatening illness or injury or permanent physical damage in connection with the use of our product or if you have to be treated as an in-patient or prolong an in-patient treatment due to the use of our product.

EN

## 3 Description

The SERVOX® Seal Pad is a sterile polyurethane foam tracheostomy compress. The underside is a light-coloured, non-adhesive layer that absorbs fluids from the tracheostoma. The SERVOX® Seal Pad can be slid over the tracheostomy tube through the perforation.

## 4 Application

- Clean the tracheostoma as usual and dry the surrounding skin carefully.
- Slide the SERVOX® Seal Pad with the lighter, fluid-absorbing side facing the patient over the cleaned tracheal cannula so that the cannula shield is completely padded. Then insert the cannula.
- To remove, remove detach the tracheal cannula from the tracheostoma, pull the SERVOX® Seal Pad downwards from the cannula and discard.

## 5 Storage, transport and disposal

Store or transport in a dry and dust-free environment at temperatures between 1°C and 30°C. Protect from direct sunlight and heat. If no disease-related instructions (e.g. MRSA, ORSA) have to be observed, dispose of the pad with household waste.

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